

SYRIS PRODUCT COMPARISON CHART

BENEFITS SEEN BY VISION ENHANCEMENT SYSTEM USERS:

- **Visualization Imm below the skin's surface**
- **Reduced reflective glare** from the skin minimizes stress, eye fatigue and potentially reduces future eye complications.
- **Increased accuracy and reduced** procedure times with vascular laser procedures, sclerotherapy, resurfacing, hair removal, hair transplantations, plastic surgery and general dermatologic procedures.
- **Reduced exam times** for pigmented lesions and other conditions.
- **Reduced recovery times** for patients (**increased accuracy** leads to reduced tissue damage from overlapping or difficulty identifying vascular structures or hair follicles).
- **Changeable magnification** on v600 and v300.
- The v600 and v300 head set models allow you to be **hands free for your procedures**.



	v600	v300	Vantage
Illumination Area	6-8 inch diffuse light	6-8 inch diffuse light	6-8 inch viewing area
Optical Barrel	<u>with</u> defined 2-3" uniform ✓ bright focal area within the center of diffuse light	<u>without</u> bright center within the diffuse area	-focal distance 8-10 in
Light Source	Halogen Light	Halogen Light	- 5 inch lens
Infrared Filter	✓ Infrared filter, lens aperture plate	No infrared filter-lens or orifice plate	- anti-reflective coating
Lasers Use	✓ laser goggles needed	If used with lasers-goggles needed	LED white lights
Comfort Bands	✓	✓	Anti reflective coating reduces ambient glare
Carrying Case	Included	available for purchase	Focal distance 8-10 inches
Photography			Pre and Post documentation capabilities
Options	Battery Pack	Battery Pack	Floor or table top photo documentation